Care management

- Robustness: excellent resistance, one-piece structure made of steel (height 80 cm) and of aluminium (height 100 cm)
- 100 % hygiene: designed with a minimum of retention zones, and resistant to disinfecting products of hospital environment
- Customization thanks to a wide range of drawers and accessories
- Ward name identification
- Ergonomic and handy: designed for medical staff comfort
- Unique design and harmonious colours









- Single-column trolley
- 2 heights: 80 cm or 100 cm
- Roller shutter door or door wing
- Code lock (with spare key) or badge lock
- 2 available inside panel models (slides) according to internal fittings







Standard configurations - PREMIUM

Structure WITH internal fittings and peripheral accessories.

Rails position:

- top right
- bottom right
- bottom left

JOINT DRAWERS - 1 small / 3 medium / 1 large

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 1 drawer 600 x 400 height 75 mm
- 3 drawers 600 x 400 height 155 mm
- 1 drawer 600 x 400 height 235 mm
- 3 standard rails
- 1 standard support for sharp container
- 1 double bag holder of 12 L
- 1 large side basket



Reference	Description
8YS64H100TJ1_RM	NEOP "Premium" 600x400 nursing trolley joint drawers H100 1S / 3M / 1L – magenta pink
8YS64H100TJ1_BC	NEOP "Premium" 600x400 nursing trolley joint drawers H100 1S / 3M / 1L – sky blue
8YS64H100TJ1_VP	NEOP "Premium" 600x400 nursing trolley joint drawers H100 1S / 3M / 1L – apple green

Label covers to order if necessary

JOINT DRAWERS – 3 small / 2 medium / 1 large

600 x 400 nursing trolley H100 equipped with:

- Sliding telescopic tray
- Door wing
- Code lock
- 3 drawers 600 x 400 height 75 mm
- 2 drawers 600 x 400 height 155 mm
- 1 drawer 600 x 400 height 235 mm
- 3 standard rails
- 1 standard support for sharp container
- 1 double bag holder of 12 L
- 1 large side basket



Reference	Description	
8YS64H100TJ2_RM	NEOP "Premium" 600x400 nursing trolley joint drawers H100 3P / 2M / 1G - magenta pink	
8YS64H100TJ2_BC	NEOP "Premium" 600x400 nursing trolley joint drawers H100 3P / 2M / 1G - sky blue	
8YS64H100TJ2_VP	NEOP "Premium" 600x400 nursing trolley joint drawers H100 3P / 2M / 1G - apple green	

Label covers to order if necessary



JOINT DRAWERS – 4 medium / 2 large

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 4 drawers 600 x 400 height 115 mm
- 2 drawers 600 x 400 height 155 mm
- 3 standard rails
- 1 standard support for sharp container
- 1 double bag holder of 12 L
- 1 large side basket



Reference	Description	
8YS64H100TJ3_RM	NEOP "Premium" 600x400 nursing trolley joint drawers H100 4M / 2G - magenta pink	
8YS64H100TJ3_BC	NEOP "Premium" 600x400 nursing trolley joint drawers H100 4M / 2G - sky blue	
8YS64H100TJ3_VP	NEOP "Premium" 600x400 nursing trolley joint drawers H100 4M / 2G - apple green	

Label covers to order if necessary

TELESCOPIC DRAWERS – 5 medium

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic top shelf
- Door wing
- Code lock
- 5 telescopic drawers F600 H160 with perforated ISO tray H100
- 3 standard rails
- 1 standard support for sharp container
- 1 double bag holder of 12 L
- 1 large side basket



Reference	Description	
8YS64H100TT1_RM	NEOP "Premium" 600x400 nursing trolley H 100 telescopic drawers 5M - magenta pink	
8YS64H100TT1_BC	NEOP "Premium" 600x400 nursing trolley H100 telescopic drawers 5M - sky blue	
8YS64H100TT1_VP	NEOP "Premium" 600x400 nursing trolley H100 telescopic drawers 5M - apple green	



TELESCOPIC DRAWERS - 1 small / 3 medium / 1 large

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 1 telescopic drawer F600 H80 with standard ISO tray H50
- 3 telescopic drawers F600 H160 with perforated ISO tray H100
- 1 telescopic drawer F600 H240 with perforated ISO tray H200
- 3 standard rails
- 1 standard support for sharp container
- 1 double bag holder of 12 L
- 1 large side basket



Reference	Description
8YS64H100TT2_RM	NEOP "Premium" 600x400 nursing trolley H100 telescopic drawers 1P / 3M / 1G – magenta pink
8YS64H100TT2_BC	NEOP "Premium" 600x400 nursing trolley H100 telescopic drawers 1P / 3M / 1G - sky blue
8YS64H100TT2_VP	NEOP "Premium" 600x400 nursing trolley H100 telescopic drawers 1P / 3M / 1G - apple green

TELESCOPIC DRAWERS - 3 small / 2 medium / 1 large

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 3 telescopic drawers F600 H80 with standard ISO tray H50
- 2 telescopic drawers F600 H160 with perforated ISO tray H100
- 1 telescopic drawer F600 H240 with perforated ISO tray H200
- 3 standard rails
- 1 standard support for sharp container
- 1 double bag holder of 12 L
- 1 large side basket



Reference	Description
8YS64H100TT3_RM	NEOP "Premium" 600x400 nursing trolley H100 telescopic drawers 3P / 2M / 1G – magenta pink
8YS64H100TT3_BC	NEOP "Premium" 600x400 nursing trolley H100 telescopic drawers 3P / 2M / 1G – sky blue
8YS64H100TT3_VP	NEOP "Premium" 600x400 nursing trolley H100 telescopic drawers 3P / 2M / 1G – apple green



Standard configurations – BASIC

Structure WITHOUT internal fittings or peripheral accessories.


TROLLEY WITHOUT FOLDABLE EXTENSION SHELF

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 3 standard rails (top right, bottom right, bottom left)

Trolley with standard inside panels – 40 mm level

Reference	Description
8YCHS164H100_RM	NEOP "Basic" 600x400 nursing trolley H100 door wing + code lock – magenta pink
8YCHS164H100_BC	NEOP "Basic" 600x400 nursing trolley H100 door wing + code lock – sky blue
8YCHS164H100_VP	NEOP "Basic" 600x400 nursing trolley H100 door wing + code lock – apple green



Trolley with inside panels for telescopic drawers

Reference	Description
8YCHST164H100_RM	NEOP "Basic" 600x400 nursing trolley H100 door wing + code lock - magenta pink
8YCHST164H100_BC	NEOP "Basic" 600x400 nursing trolley H100 door wing + code lock - sky blue
8YCHST164H100_VP	NEOP "Basic" 600x400 nursing trolley H100 door wing + code lock - apple green

TROLLEY WITH FOLDABLE EXTENSION SHELF

600 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 2 standard rails (bottom right, bottom left)
- 1 foldable extension shelf

Trolley with standard inside panels – 40 mm level



Reference	Description
8YCHS164H100T_RM	NEOP "Basic" 600x400 nursing trolley H100 H100 door wing + code lock - foldable extension shelf – magenta pink
8YCHS164H100T_BC	NEOP "Basic" 600x400 nursing trolley H100 door wing + code lock - foldable extension shelf — sky blue
8YCHS164H100T_VP	NEOP "Basic" 600x400 nursing trolley H100 door wing + code lock - foldable extension shelf – apple green

Trolley with inside panels for telescopic drawers

Reference	Description
8YCHST164H100T_RM	NEOP "Basic" 600x400 nursing trolley H100 door wing + code lock - foldable extension shelf – magenta pink
8YCHST164H100T_BC	NEOP "Basic" 600x400 nursing trolley H100 door wing + code lock - foldable extension shelf – sky blue
8YCHST164H100T_VP	NEOP "Basic" 600x400 nursing trolley H100 door wing + code lock - foldable extension shelf – apple green



Custom configurations

STEP 1 - Trolley structure

Choice of trolley structure according to inside panel model

	Trolley Height 80 cm*	Trolley Height 100 cm*
Standard inside panels 40 mm level	8Y164H80P4	8Y164H100P4
Inside panels for telescopic drawers	-	8Y164H100T


This reference includes:

- Structure
- Wheels
- Inside panels

STEP 2 - Trolley customization Selection of elements 1 to 4 is compulsory		-	Trolley Height 80 cm*	Trolley Height 100 cm*
		thermoplastic	8Y001_6400TH	8Y001_6400TH
1	Upper shelf	compact laminate	8Y001_6400CP1	8Y001_6400CP1
		Roller shutter door	8Y002_6008F3	8Y002_6010F3
2	Closing system	Door wing for standard inside panels	8Y002_0008F1	8Y002_0010F1
	Зузст	Door wing for telescopic drawers	-	8Y002_0010F2
		For roller shutter door	8Y003CODE3	8Y003CODE3
	Code lock	For door wing for standard inside panels	8Y003CODE1	8Y003CODE1
3		For door wing for telescopic drawers	-	8Y003CODE2
•		For roller shutter door	8Y003BADGE3	8Y003BADGE3
	Badge lock	For door wing for standard inside panels	8Y003BADGE1	8Y003BADGE1
		For door wing for telescopic drawers	-	8Y003BADGE2
		Magenta pink	8Y004_0008RM	8Y004_0010RM
4	Colour	Apple green	8Y004_0008VP	8Y004_0010VP
		Sky blue	8Y004_0008BC	8Y004_0010BC

^{*}For precise height of working surface, see technical data







Roller shutter closure



Door wing



Code lock



Badge lock



Accessories and internal fittings

Extension shelf, standard rail, ergonomic handle, motorization, arm for laptop, upper frame, drawer... *Please refer to product data sheets.*















Technical data

Trolley with thermoformed plastic top shelf

Cart with top shelf in compact laminate





The overall dimensions are measured with the wheels in the least bulky position. The depth incorporates a lock by default.





Materials

Structure	Complete welded one-piece structure, EPOXY anti-corrosion treated H 80: steel (electro-galvanized) / H 100: aluminium		
Work surface	Thermoplastic: ABS / PMMA or compact laminate 10 mm thickness		
Slides	ABS (Acrylonitrile Butadiene Styrene) or ballbearing slides made of galvanized steel with shock absorber		
Closing system	Roller shutter door: PP (Polypropylene), 2 mm thickness Door wing: Aluminium with EPOXY painting		
Wheels	Simple non marking rubber wheels - Polyamide cover (PA)		
Bumpers	4 non marking round stop blocks made of EVA (Ethylene Vinyl Acetate)		
Strip of colour	70 micron calendered PVC vinyl film, with acrylic adhesive		

Additional information

Wheels	125 mm diameter - 2 breaking wheels + 2 swivelling wheels - wheel and wheel-pin double blocking		
Empty weight	H 80: 27 kg	H 100: 24 kg	
Capacity	H 80: - 15 levels of 40 mm for trolley with roller shutter door - 16 levels of 40 mm for trolley without roller shutter door or trolley with door wing	H 100: - 19 levels of 40 mm for trolley with roller shutter door - 20 levels of 40 mm for trolley without roller shutter door or trolley with door wing - 10 telescopic drawers	
Locking system	Code lock / Badge lock		

Packaging

Packing	Individually on pallet strapped and film-	Dimensions	H 80: 80 x 120 x 120 cm
	wrapped, or on wheels, with proper protection	L X D x H	H 100: 80 x 120 x 145 cm
Gross weight (trolley + packaging)	Н 80: 39 kg Н 100: 36 kg	Volume	1,15 m ³



