











# **BLUE TRANSFER BOARD**

# The safe way to transfer patients

A legal requirement in some countries, taken for granted in others – working conditions should be as comfortable as possible. Healthcare workers should avoid lifting when transferring patients. SAMARIT Products assist with transfer, without the need to lift. The products are designed to facilitate the transfer without using a lot of effort, thus making the work easier for the carers and kinder to their backs. Can be ordered as OEM product.

**BLUE LONG WIDE BLUE SHORT NARROW BLUE LONG NARROW BLUE SHORT WIDE** 

EU-SRN: CH-MF-000016280 CHRN: CHRN-MF-20000184







#### **Samarit Medical AG**

Gewerbestrasse 12 CH-8132 Egg, Schweiz

Tel.: +41 44 918 10 11 Fax: +41 44 918 24 39 e-mail: info@samarit.com - www.samarit.com





## The BLUE TRANSFER BOARD

#### **BLUE LONG WIDE**

foldable



Suitable for radiology, operating theatre, intensive care units and all special wards.

**Size**: 180 x 50 cm / 89 x 50 cm, 2.9 kg 70 x 20 inch / 35 x 20 inch 6.4 lbs

SWL\*: 240 kg / 529 lbs

Order No. 440.8300 UDI/GTIN: 7640144460228 **BLUE SHORT NARROW** 

non-foldable



For ambulance, operating theatres, wake-up stations and emergencies.

**Size**: 100 x 42 cm, 1.4 kg 39 x 17 inch 3.1 lbs

**SWL\***: 150 kg / 331 lbs

Order No. 440.8301 UDI/GTIN: 7640144460235

#### **BLUE LONG NARROW**

foldable



Suitable for operating theatre, X-Ray, various wards.

**Size**: 180 x 42 cm / 89 x 42 cm, 2.1 kg 70 x 17 inch / 35 x 17 inch 4.6 lbs

SWL\*: 180 kg / 397 lbs

Order No. 440.8305 UDI/GTIN: 7640144460259

\*SWL = Safe Weight Limit

# 4 different models to suit every hospital application

- High performance
- Radioluscent
- Non-toxic
- Anti-bacterially surface
- Safe and comfortable for patients
- Prevents back injuries to care workers
- Skin-friendly (no latex)
- Easy to clean and disinfect hygiene approved
- Disposable covers available
- Tear resistant nylon cloth
- Bridges gaps up to 20 cm
- Lightweight, easy suspendable

### **BLUE SHORT WIDE**

non-foldable



For urology operating theatres, gynaecology operating theatres, wake-up stations and emergencies.

**Size**: 80 x 50 cm, 1.4 kg 3.15 x 2 inch 3.1 lbs

SWL\*: 150 kg / 331 lbs

Order No. 440.8304 UDI/GTIN: 7640144460242