

NURSING HOME
HOMECARE
AMBULANCE



ECONOMIC ROLLBOARD H-LINE

The safe way to transfer patients

Healthcare workers should avoid lifting when transferring patients. The products are designed to facilitate the transfer without using a lot of effort, thus making the work easier for the carers. Can be ordered as OEM product.

H-LINE 170
H-LINE 120

The ECONOMIC ROLLBOARD H-LINE

H-LINE 170

foldable



Medium use for Nursing Homes, Homecare and Ambulances.

Size: 170 x 45 cm / 83 x 45 cm, 2.4 kg
67 x 18 inch / 33 x 18 inch

Material:
Sliding fabric around a reinforced polyethylene core.

Order No. 440.3170

H-LINE 120

foldable



Medium use for Nursing Homes, Homecare and Ambulances.

Size: 120 x 45 cm / 60 x 45 cm 1.3 kg
47 x 18 inch / 24 x 18 inch

Material:
Sliding fabric around a reinforced polyethylene core.

Order No. 440.3120

2 different models to suit care and ambulance

- Good performance
- Radioluscent
- Anti-bacterial PU surface
- Comfortable for patients
- Skin-friendly (no latex)
- Easy to clean and disinfect
- Washable cover 95°
- Disposable covers available
- Lightweight
- 1 year warranty